



Warming or cooling a patient can be a critical factor in achieving positive patient outcomes. Rapr•Round™ Body Wraps leverage existing circulating water technology combined with a unique design that literally wraps around the patient. Compared to conventional water blankets this design increases the area of body surface contact and, therefore, is a more efficient warming/cooling system.

The Rapr•Round is effective for inducing and maintaining mild hypothermia, rewarming and temperature management of febrile patients. Used with the Gaymar Medi-Therm® hyper/hypothermia unit, the system provides clinically effective therapy at a cost effective price.

RAPR•ROUND™

Hyper/Hypothermia Body Wraps

DHV530 – Small/Medium Vest DHV535 – Large Vest DHL540 – Leg

Rapr•Round™ for the management of hyper/hypothermia

Features

Up to 72% of the body is in contact with the blanket

Conductive (water) Technology

Non-woven fabric against the patient

Two sizes of Vests

One size fits all Legs

Vest separate from the Legs

Easy to apply

Benefits

- Greater effectiveness in warming and cooling
- Increased skin surface contact, greater than standard modalities
- Increased heat transfer

- As a heat conductor, water is four times more effective than air

- Gentle on the skin
- No sheet between the Body Wrap and the patient required

- Will cover up to 90% of the adult patient population
- Only two items to stock

- Ease of use, only necessary to stock one size of Leg Wrap, no right or left

- Allows for flexibility in using the product
- Allows for easy patient access
- One product can be used in multiple hospital departments

- No adhesives to adhere to the skin
- Simple Velcro® style closure
- Wraps move with the patient when repositioning occurs

Ordering Information

Catalog Number	Description	Packaging
DHV530	Small / Medium Vest chest sizes 32–46 inches	5 Vests per case
DHV535	Large Vest chest sizes 46 inches and greater	5 Vests per case
DHL540	Leg (right or left) one size fits all (for adults)	10 Legs (5 sets) per case

Single patient use • Latex free

For a complete patient warming and cooling system, use Rapr•Round Body Wraps with the Gaymar Medi-Therm® hyper/hypothermia unit.



Gaymar provides a comprehensive range of Temperature Management and Pressure Ulcer Management products. To order or inquire, call Gaymar's Customer Service Department at 1.800.828.7341.

Gaymar Industries, Inc.

10 Centre Drive

Orchard Park, NY 14127

Toll Free 1.800.828.7341

Fax 1.800.993.7890

International +1 716.662.8636

International Fax +1 716.662.0730

www.gaymar.com

Available in Canada through

Waterloo Bedding Co. Ltd.

825 Trillium Drive

Kitchener, ON N2R 1J9

Toll Free: 1.800.203.4293 • 519.742.4447

Fax: 519.742.3699

Gaymar and Medi-Therm are registered trademarks of Gaymar Industries, Inc.

Velcro is a registered trademark of Velcro Industries B.V.

©2006. Gaymar Industries, Inc. MAW4928 R8-06

 **GAYMAR®**
Temperature Management